


Light & Ergonomic Pipette

# NICHIRYO *Le*

**New NICHIRYO Pipette**  
**Lightens your stress during**  
**long hours of daily work.**  
**Surprisingly light in weight and**  
**cost effective.**



- ◆ Only 83g in weight (NLE-100)
- ◆ User-friendly ergonomic design
- ◆ Excellent cost performance

BE ORIGINAL

# NICHIRYO

# Light, Easy-to-use and Cost effective. This is completely new pipette from Nichiryo.

## Feel differences of **NICHIRYO Le**

### ◆The grip is different

Enlarged and stream lined grip fits well with your fingers,  
That can be kept in hand firmly and be used with much less stress.

### ◆The push-button is different

The user-friendly enlarged push-button can be pressed easily and can  
be used for long time with less probability of RSI(repetitive strain injury).

### ◆The weight is different

The weight is just 83g(NLE-100).This weight is perfect for a  
tool of the daily work and reduce probability of RSI for  
long time operation.

### ◆The volume setting is different

No need to lock the dial after volume setting for your immediate  
operation.Easy to set a volume with the enlarged push button.

### ◆The ejector is different

No screw is used for the tip ejector.  
The tip ejector is one-touch detachable and attachable.

### ◆The price is different

Completely new designed and drastic cost-down pipette with  
unchanged Nichiryo quality.

### ◆The after-sales service is different

NICHIRYO Le is not disposable pipette.  
Continued after-sales service is given as promised as other  
Nichiryo pipette.

## » NICHIRYO Le Line Up



The pipetting volume covers from 0.5 $\mu$ L to 1000 $\mu$ L,  
Widely used volume range for research works.

Cat.No.		NLE-10	NLE-20	NLE-100	NLE-200	NLE-1000
Volume range		0.5~10 $\mu$ L	2~20 $\mu$ L	10~100 $\mu$ L	20~200 $\mu$ L	100~1000 $\mu$ L
AC%	Measured (max.)	$\pm 1.2$	$\pm 1.2$	$\pm 1.0$	$\pm 0.8$	$\pm 0.7$
	Measured (min.)	*	* $\pm 5.0$	$\pm 2.0$	$\pm 1.2$	$\pm 1.2$
CV%	Measured (max.)	$\leq 0.8$	$\leq 0.8$	$\leq 0.6$	$\leq 0.6$	$\leq 0.5$
	Measured (min.)	*	* $\leq 3.0$	$\leq 1.0$	$\leq 1.0$	$\leq 0.8$
Weight		83.4g	83.8g	83.0g	83.3g	83.3g
Usable tips	Bulk	BMT-SS	BMT-SG	BMT-SG	BMT-SG	BMT-L
	Racked	BMT-SSR	BMT-SGR	BMT-SGR	BMT-SGR	BMT-LR

\* The AC and CV values at 0.5 $\mu$ L by NLE-10 and 2 $\mu$ L by NLE-20 largely depend on the operator's skill and the environment under which the pipette is used.

\* The AC and CV values are the value by use of the disposable tip described in this leaflet.

\* The specification and price are subjected to change due to modification for quality improvement without notice.

<http://www.nichiryo.co.jp>

Distributor :

 **NICHIRYO CO.,LTD.**

4-10, IWAMOTO-CHO 2-CHOME, CHIYODA-KU, TOKYO 101-0032, JAPAN  
TEL.81-3-5820-2664 FAX.81-3-5820-2410

E-mail [i-sales@nichiryo.co.jp](mailto:i-sales@nichiryo.co.jp)